



Main System							for stockings / bandages	
Models	1ch/Handy	2ch/Extend	5ch/Desktop	10ch/Desktop	1ch/Desktop	1ch/Desktop and many points are possible		
	AMI 3037-2	AMI 3037-2B	AMI 3037-5S-II	AMI3037-10-II	A0101-G35k-II	A0905-SA-II	AMI3037-SB-nP-II	AMI3037-SB-mH-II
Unit of pressure	convert the voltage	convert the voltage	k Pa		k Pa	k Pa (or mmHg)	h Pa	(mmHg)
ch	1	2	5	10	1	1	1	1
Measuring Range	0~34.00 kPa		0.00~35.00 kPa		0.0~35.0kPa	0.0~35.0 kPa	0~350 hPa	0~263 mmHg
Air-pack (φ20)	—		—		—	0.0~20.0 kPa	0~200 hPa	0~150 mmHg
Measuring Accuracy in 23°C Environment	0~14.00/±0.1 kPa 14.00~34.00/±0.25 kPa		0~14.00 ±0.2 kPa 14.00~35.00 ±0.3 kPa		0~14.00 ±0.2 kPa 14.00~35.00 ±0.3 kPa		±2 hPa	±2 mmHg
Output	1 kPa → 0.1V (100mV)				1 kPa → 0.01V(10mV)		—	—
Numerical display	—	—	LED 5ch	LED 10ch	LED 1 ch	LED 1 ch	LED 1 ch	LED 1 ch
Durable Pressure	up to 70 kPa						up to 700 hPa	up to 500 mmHg
Power supply	AC100V(11VA) AA battery/1.5v×4	supplied by the left instrument	AC90~250V 22.5VA 45VA		AC100V 30VA	AC100V 30Va	AC100V 30VA	
Air-Cylinder	Optional attachment : AMI3037-AC/ Black, red, green, and blue					Equipped by machinery	Equipped by machinery	
Air-pack Which is connected	Standard:AMI3037-P2,AMI3037-P5/Parasol type : AMI3037-P2P,AMI3037-P5P					φ 20/?m AMI3037-SB-P3/P5	AMI3037-SB-P3/P5(φ 20/?m)	
Zero adjustment	Adjustment by screw		Reset switc/hzero switch		Automatically		Automatically	
Applicable Environment	0~45°C RH. 25~85%						0~45°C RH. 25~85%	
External Dimensions	W92×H45×D160		W510×H240×D230		W245×H65×D175		W245×H65×D175	
Weight	0.6Kg	0.6Kg	8Kg	10Kg	2Kg		2Kg	
Feature	Small and a Battery	Small and a Battery	Measure 5ch at the same time	Measure 10ch at the same time	The consecutive output	The consecutive output/Many points are possible	Watch numerical value/Many points are possible	

Air-Pack Model-Standard: AMI3037-P2 (Two) ,AMI3037-P5 (Five) /Parasol type: AMI3037-P5P (Five) ,AMI3037-P10P (Ten)									
Product number	①	②stock size	③	④	⑤	⑥	⑦Finger use	⑧⑨/⑩	
Size (mm)	φ 15	φ 20	φ 25	φ 30	φ 8×L=28	φ 20×L=30	φ 12×L=15, φ 15×L=17	⑧⑨Tube is connected under a bag or set in the bag.	
Size & Tube Length	φ15/1.5m-2m	φ 20/2m-3m	φ 25/2m-3m	φ 30/2m-3m	I /1.5m-2m	G/2m-3m	φ 12×15/1.5m		
Parasol type(φOP)	—	φ 20P/2.5m	φ 25P/2.5m	φ 30P/2.5m	T /1.5m-2m	GP/2.5m	φ 15×17/1.5m		
for A0905 and AMI3037-SB	—	SB φ 20/1m - 0.5m	Model A0905: Can use the special air-pack						—
Max. Measuring range (Approx.)	10	20	30	38	15	20	15		
Measuring accuracy (sensitivity to soft material)	±0.2	±0.1	±0.1	±0.1	±0.25	±0.1	±0.3	—	
Size Cover tape/Models	φ 47/AMI3037-PTS		φ 47/AMI3037-PTB		—				
Air Cylinder	Black(half of Red)	Red	Green	Blue	Black	Red	Black	—	
Main Applications (mm)	Sphere of φ 15 or larger	Sphere of φ 20 or larger (standard)	Sphere of φ 25 or larger but with small curvature	Sphere of φ 30 or larger but with small curvature	Objects of 8 in width and 28 in length or longer. For belts, flat cords, side of fingers.	Objects of 20 in width and 30 in length or longer. For averaged portion of shoulder of coats.	Objects of 12 in width and 15 in length or longer. For fingers, etc.	—	



Pressure Calibration Unit	
Models	AMI3037-CS
Calibration Range	0~4.5 kPa
Compositions	
①Fine Chemical Cylinder	φ 95 × L=500
②Scale	L=500
③T-shaped bar	φ 50 L=500
④Pipe with tube	φ 6 × L=500 1.5m
Weight	1kg



For the abnormal management of the measurement



AMI Techno co.,ltd.

fax: +81-3-5339-7414 3-5-3-1313, NISHI-SHINJUKU,
SHINJUKU-KU, TOKYO, 160-0023, JAPAN
E-mail: ami-tec@m2.pbc.ne.jp <http://www.ami-tec.co.jp>

« THE EXCLUSIVE OVERSEAS DISTRIBUTOR »


SANKO TSUSHO co.,ltd.

tel: +81-3-5777-3627 fax: +81-3-5777-3629
TOKYO, JAPAN E-mail: sales@sankotsusho.co.jp



Air-pack type contact surface pressure continuously measuring system

Air-pack type Air packed in an air-bag is put between two contacting materials and air-pressure brought about inside that air-bag is measured.



Measurement of body pressure and of pressure in its sitting attitude.

Measurement of pressure appearing in clothes.




Usage

Contact pressure



Restraint pressure

Body pressure

Clothes pressure



The pressure taking place on the hipbones will be 20kPa with a hard plane. How may be the conditions or curve or cushion?



New

Living body (soft condition) can be measured.

You can measure just by pulling the knob

A0905-SA-II Both the living body (soft condition) and the dummy (hard condition) can be measured.

Pressure by sitting

Flexibility

The accuracy of the measured value can be certified by the water pressure calibration method.

Surface Pressure
Blood Flow
Skin Temperature




Patented

U.S.A.: U.S.Pat.5522966 July 13,1999
Europe: E.P.Pat.0831313 August 9,2000
JAPAN 3268601 January 18,2002

Air enclosure method contact pressure/PT#,5601489
Cylinder-shaped expansion-style contact pressure calibration method/PT#,5522337
Pressure sensor, built-in system /PT#,5601489

... Measured by water pressure calibration method ...



Calibration

Molding of foaming materials

«MANUFACTURER»

AMI Techno CO.,LTD. TOKYO,JAPAN

«THE EXCLUSIVE OVERSEAS DISTRIBUTOR»

SANKO TSUSHO CO.,LTD. TOKYO,JAPAN

Tel: +81-3-5777-3627 E-mail: sales@sankotsusho.co.jp


This contact surface pressure measuring system is a system to measure a surface pressure of a soft material to which the sensor of this system is in contact and to correct errors in the reading to be brought about by a thermal change in the environment. This system has been approved by patent authorities in several countries in the world as a unique system to be able to effect a reliable measurement of a contact pressure of soft materials with a high reproducibility. This system does not require any special preliminary correction for the reading to be usually made before use. This system can also make a continual measurement sequentially.

Air-pack method

An air is enclosed between two contacted surfaces and the pressure of that enclosed air is measured externally from the contacted surface. Into a flat bag made with a flexible film of a least possible elasticity, an air is enclosed in the thickness of 1mm.

Feature

- (1) It is easy to stick each part of the living body and between ruggedness of the cloth and reproducibility is good in 20mm (standard) diameter. The center part of the bag comes in contact with the point when it is too large and the situation is changeable and reproducibility worsens partial pressure it. For example, only the center of the air pack comes in contact in the ankle that the ruggedness is 30mm in diameter.
- (2) The influence of the air pack that infuse some air in 1mm in thickness becomes thinner because of pressure sinks softly of the living body and the clothes is extremely little. However, when infusing in 1.5mm in thickness and measuring it in the cylinder of $\phi 100$, the value about 1.3 times 1mm in thickness was confirmed.
- (3) The circle shape mounted to body part easily. The result change by the direction when designed besides the circle. For example, the top becomes the mean value of the lowering area high in surroundings in the measurement with the hemisphere.
- (4) In piping in the tube of 1mm that doesn't collapse, it is unaffected in the piping situation catching the transformation of the bag of the air pack (alteration in volume). The result change for a soft material such as silicon rubber in the piping situation.



Each error in measurement




①Permissible error+②Measuring area error+③Position error+④Restoration of tools material+⑤Attitude error

- ①**Permissible error:** Describes in the specification.
- ②**Measuring area error:** The difference of the characteristic the mounted position, the figure and the softness of the air pack sensor is caused.
- ③**Position error:** The value changes partially in how like hitching and a horizontal gap to match it when setting it to the measurement site when tools are installed. The error is somewhat caused though putting the sign and the line in tools in detail, and improving accuracy are necessary.
- ④**Restoration of tools material:** The difference of the restoration is caused by the sweat, the expansion in the rig frequency.
- ⑤**Attitude error:** It influences pressure in transmission and the weight shift of the material tension because of the difference among the posture angle, arm and the head position.

Note: In the measurement of body pressure to know the distribution pressure is important. However, the method of paving the seat sensor in the distribution pressure measurement and the calling contact surface and the measurements has the anxiety in reproducibility and the measurement accuracy. Non-elastic film seat enclosed between two contacted surfaces and A flexible characteristic and the contact friction are completely changed and it is not easy to call body

Calibration for contact pressure and heat influence

The plane pressure can be calibrate by mounting the air pack to horizontal plate. The contact surface is hardness for horizontal plate and softness (living body) for water side.




Water thermally controlled by a water heating apparatus is poured into a vessel and is held at a constant water level. The Air-Pack is fixed horizontally in water in the depth of 400 mm.

[Result] change 40°C is shown to be 0.1 kPa

<The reasons why a thermal effect is so small are as follows:>
 ① Although a thermal expansion factor of air is 1/(273+atmospheric temperature), a voluminal swell of the Air-Pack to thickness direction is only a little because its shape is flat.
 ② Because an amount of the enclosed air is smaller than the maximum voluminal capacity of the Air-pack, there is a sufficient remaining capacity inside the Air-pack to be able to absorb a swell of the enclosed air.
 ③ In actual measurement, a swell to the thickness direction of the

Calibration for pressure sensor and output

A pressure on horizontal plane corresponding to the water pressure can be got by measuring the water depth.



SI unit		Former unit			Model	AMI3037-2, 5S, 10		A0101 / A0905	
hPa	kPa	(gf/cm ²)	(mmH ₂ O)	(mmHg)		DC output		DC output	
						DC V	DCmV	DC V	DCmV
10	1	10.2	102	7.5		0.100	100	0.010	10
0.98	0.098	1	10	0.736		0.010	9.8	0.001	0.98
0.098	0.0098	0.1	1	0.074		0.001	0.98	0.0001	0.1
1.333	0.1333	1.36	13.6	1		0.013	13.3	0.0013	1.3



Purpose of use for pressure measuring device

The measurement accuracy and reproducibility are confirmed the reproduction of a flexible plane contact environment. Error margin confirmation of secular distortion and confirmation of unexpected value when measuring it. It is not necessary to calibration before each measurement. * Because the standard of the contact pressure (ISO and JIS, etc.) in soft respect is not decided, the proofreading examination such as the inspecting agencies cannot be done under the present situation.

System

An air infuse to the air-pack sensor with air cylinder. And it connect to main unit and air-pack sensor mount to the body.

The air-pack sensor connect to the body in advance and wear the clothes. An air infuse to air-pack sensor with internal air cylinder by pulled the cock.



(former) Accessory: Air cylinder

2 type measurement

● <1ch·3ch·5ch·7ch...> Handy type/Extension AMI3037-2 (1ch) / Extension unit AMI3037-2B (2ch)


Main unit internal model/Semi-auto air cylinder

Feature: 2 type power supply / The output signal is DC voltage and small main unit with dry battery (4 AA batteries) and AC adaptor (AC100v). It automatically changes to the AC adaptor when the battery consumed.

Add channel/ Main unit AMI3037-2 is 1 ch and it can add each more 2 ch / 1 unit with extension unit AMI3037-2B.

Composition: ① Main unit/AMI3037-2(1ch), ② Extension unit AMI3037-2B(2ch), ③ Air cylinder/AMI3037-CY ④ Air-pack sensor 5pcs/AMI3037-P5, ⑤ Covering tape(200sheet)/AMI3037-PT, ⑨ AC adaptor/AMI3037-AC ● Calibration set/AMI3037-CS accessory: ⑥ Press tester ● screw driver for zero

Option: ● Data logging device fo PC AM8051 (6ch) Accessory : software ● USB communication cable ● A C adapter (for data logger)



● <5ch·10ch> Main unit AMI3037-5S-II (5ch) AMI3037-10-II (10ch)

Accessory: Air cylinder



Feature: Main unit (10 ch desktop type): For multi-point measurement with DC voltage output, it can measure the multi-point at one time and the data transfer to the data logger.

Compositio ① Main unit/AMI3037-10(10ch) ③ Air cylinder/AMI3037-AMI3037-CY ④ Air-pack sensor n : 5pcs X 3set/AMI3037-P5 X 3 ● Covering tape(200sheet)/AMI3037-PT ● Calibration set/AMI3037-CS, accessory: ⑥ Press tester ⑦ screw driver for zero adjust

Option : ● Data logging device for PC (10ch or more)

● <1ch> Main unit with display (Auto zero adjust) A0101-G35-AC (1ch)

Accessory: Air cylinder

Feature: Auto Zero Adjust: 1 ch with the output signal (DC voltage output).Auto zero adjustments after turned on the power. AC100V power supply.


Composition: ① Main unit/A0101-G35-AC(1ch) ④ Air-pac sensor 3pcs/SB-φ 20/0.5m·1m ⑤ covering tape(200sheet)/AMI3037-PT ⑩ output cable ● Calibration set/AMI3037-CS, accessory: ⑥ press tester(small)

Option: ● Data logging device fo PC AM8051(6ch) Accessory : software ● USB communication cable ● AC adapter (for data logger)



● <1 ch> Main unit with display /DC output/Auto zero adjust) A0905-SA-35k (1ch)

semi-auto air cylinder



For multipoint : Air-pack sensor previously mounted to the part and measurement tools install to the body. And the air infuse by the built-in air cylinder. Only the number of air packs repeats the measurement. Moreover, the change can be measured by being reconnected and measuring it after time passes with the air pack mounted.

Semi-auto air cylinder 2type<Black/Red> : It can select that air amount of infusion enclosed by the shape of air-pack and hardness on the measurement side. I case of Air-pack sensor φ 20, pull the red cock for soft living body with thickness 1mm and pull the black cock for hard materials with thickness 0.5mm.


Auto Zero Adjust : 1 ch with the output signal (DC voltage output).Auto zero adjustments after turned on the power.

Composition: ① Main unit/A0905-SA-35k(1ch) ⑨ relay tube 1. 5m ④ Air-pac sensor 3pcs/SB-φ 20/0.5m·1m ⑤ covering tape(200sheet)/AMI3037-PT ⑩ output cable ● Calibration set/AMI3037-CS, accessory : ⑥ press tester(small)

Option: ● Data logging device fo PC AM8051(6ch) Accessory : software ● USB communication cable ● AC adapter (for data logger)

● <1ch> Mai unit with display (for multipoint) AMI3037-SB-hP (hPa) /mH (mmHg) (1 ch) <stocking/bandage>

semi-auto air cylinder



Feature: For multipoint : Air-pack sensor previously mounted to the part and measurement tools install to the body. And the air infuse by the built-in air cylinder. Only the number of air packs repeats the measurement. Moreover, the change can be measured by being reconnected and measuring it after time passes with the air pack mounted.

Auto Zero Adjust : 1 ch with the output signal (DC voltage output).Auto zero adjustments after turned on the power.

① Main unit/AMI3037-SB(1ch), ⑨ relay tube 1. 5m ④ Air-pack sensor 5pcs/SB-φ 20/0.5m·1m, ⑤ Covering tape(200sheet)/AMI3037-PT, accessories ⑥ press tester(small)